


		Minitom
Capacity	Height x Length	30 x 40 mm (1.2" x 1.6")
	Diameter	40 mm (1.6")
	Maximum cutting force	0.34 N
Cut-off wheel	Diameter	100 - 127 mm (4" - 5")
	Arbor diameter	12.7 mm (0.5")
Cut-off wheel motor	Rotational speed	100 - 420 rpm
Coolant tank	Volume	250 ml (8.5 fl oz)
Software and electronics	Controls	Touch pad
	Display	N/A
Safety standards		See the Declaration of Conformity
REACH		For information about REACH, contact your local Struers office. For information about REACH, contact your local Struers office.
Operating environment	Surrounding temperature	5 - 40 °C (41 - 131 °F)
	Humidity	< 95 % RH non-condensing
	Noise level	L _{pA} = 58 dB
Power supply	Voltage/frequency	100 - 240 AVC 50 - 60 Hz
	Power inlet	1L + (N) + PE
	Power S1	N/A
	Current, nominal load	1.4 A
	Current, max.	2.8 A
	Ampere rating of the largest motor or load	3.2 A
Residual Current Circuit Breaker (RCCB)		Type A, 30 mA (or better) is recommended
Noise level	A-weighted sound emission pressure level at workstations	L _{pA} = 58 dB
		Uncertainty K = 4 dB
		Measurements made in accordance with EN ISO 11202
Vibration level		N/A



TECHNICAL DATA

Minitom

Minitom		
Dimensions and weight	Width	27.5 cm (10.8")
	Depth (main body)	28 cm (11")
	Depth (with handle)	41 cm (16")
	Height (main body)	20.5 cm (8.1")
	Height (with handle)	25.6 cm (10.1")
	Weight	12 kg (26.5 lbs)

Doc. no.: 14437750_C_en

Date of release: 2026.07.30